

Product Compatibility Chart

Discipline	Market/Application	Endoscope Series	Processors / Lightsources								AquaTYPHOON / PlasmaTYPHOON / PlasmaTYPHOON+ and PlasmaBAG Systems
			EPK-i8020c INSPIRA™	EPK-i7010 OPTIVISTA plus	EPK-3000 DEFINA	EPK-i5500c IMAGINA	EPK-V1500c VERSA	CP-1000 VIVIDEO Video Processor	KAY9372HD Digital Capture	BS-LL1 LED Light Source Battery-type	
Upper GI and Lower GI	Advanced diagnosis and treatment	i20c	X								X
		i10	X	X							X
		K10		X	X						X
		90i	X	X							X
		90K	X*	X	X						X
	Diagnostic and standard therapeutics	i10c	X			X					X
	Basic diagnosis and treatment	V10c					X				X
Biliopancreatic	Advanced diagnosis and treatment	i10, i10c, J10U	X	X							X
	Diagnostic and standard therapeutics	i10c, J10U	X	X	X (J10U)						X
	Basic diagnosis and treatment	i10c, J10U	X	X	X (J10U)						X
Pulmonology	Advanced diagnosis and treatment	i10c, J10, J10U	X	X	X (J10, J10U)						X
	Diagnostic and standard therapeutics	J10, J10U, 70K	X (J10, J10U)	X	X						X
	Basic diagnosis and treatment	V, RBS								X (RBS-series)	X
Intubation	Advanced diagnosis and treatment	70K			X						X
	Basic diagnosis and treatment	RBS								X	X
ENT	Advanced diagnosis and treatment	J10		X	X						X
	Diagnostic and standard therapeutics	J10, CP		J10	J10			CP			X
	Basic diagnosis and treatment	RP3, RBS			X					RBS	X
Urology	Diagnostic and standard therapeutics	75K									X
	Basic diagnosis and treatment	RBS								X	X

*Limited compatibility