



Endopress 442

Mechanical
lymph drainage for
professional medical use





Endopress 442

The EndoPress 442 is the near-silent pressotherapy device by Enraf-Nonius to provide effective and comfortable lymph drainage treatment. Equipped with various programmes that offer variation in the pressure settings in the various segments. The built-in clinical programmes offer the most suitable settings for the various indications. This makes it possible to apply different lymph drainage techniques. The pressure can be adjusted quickly and accurately at a value between 20 and 150 mmHg for the entire cuff or for each individual segment.



Endopress 442


- The pressure can be adjusted per segment
- The double arm cuff consists of 5 segments, the single arm cuff consists of 7 segments and the leg cuff consists of 9 segments. The pressure build-up can be configured per segment.
- A very quiet and powerful compressor of the highest quality
- Easy to operate via the full colour touch screen
- Access to 21 clinical protocols
- Memory for storage of more than 100 favourite programmes
- Very comfortable for therapist and patient
- Easy to clean
- Optionally available: handy storage bag that attaches to the treatment couch for neat storage of the various cuffs.



Endopress 442

Reliable, durable and quiet

Lymph drainage – also referred to as pressotherapy – is the stimulation of the lymphatic system by means of a specific, subtle massage, allowing for the drainage of fluid and waste products that are present in and around the cells.



Endopress 442

1442921 Endopress 442 230V

1442922 Endopress 442 115V

Device for mechanical lymph drainage

ORDERING DATA



Enraf-Nonius


Endopress 442

TECHNICAL SPECIFICATIONS

Mains supply	230 VAC 50-60 Hz (1442921) 115 VAC 50-60 Hz (1442922)
Pressure settings	20-150 mbar
Programmable treatment time	60 minutes
Output connector	1
Pre-programmed treatment suggestions	21
Free programmable treatment suggestions	>100
Dimensions (HxWxD)	48 x 56 x 110 cm

This is a medical device.

It complies with all applicable medical regulations. For more details please consult the instructions for use (downloadable via www.enraf-nonus.com).



Endopress 442

ACCESSORIES

Optional accessories

- 1442802 Cuff for abdomen, legs and feet
- 1442803 Cuff for 2 arms
- 1442804 Cuff for 1 arm
- 1442808 Storage bag, with 1 zippered storage compartment on each side of the treatment couch



1442802



1442803



1442804



1442808





Endopress 442

Mechanical
lymph drainage for
professional medical use

Enraf-Nonius B.V. | Vareseweg 127 | 3047 AT Rotterdam | The Netherlands
www.enraf-nonius.com | info@enraf-nonius.nl | + 31 - (0) 10 20 30 600



PARTNER FOR LIFE