



## Portable Massage Couch

The safe,  
easily portable couch

# Portable Massage Couch

A practical and safe alternative for applying a treatment, wherever you are. Easy to transport and quick to set up.





## Portable Massage Couch


- Weighs only 15 kg
- Load capacity up to 135 kg
- Attachable head section with nose hole
- 56 cm wide sections
- Aluminium frame with foldable legs
- Fixed height of 73 cm



# Portable Massage Couch

3450007 Portable couch with headsection


## ORDERING DATA



## Portable Massage Couch

Size tabletop (lxw):	210x56 cm
Fixed height:	73 cm
Load capacity:	135 kg
Upholstery:	grey
Color frame:	aluminum
Weight:	15 kg

## TECHNICAL SPECIFICATIONS





# Portable Massage Couch

The safe,  
easily portable couch

Enraf-Nonius B.V. | Vareseweg 127 | 3047 AT Rotterdam | The Netherlands  
[www.enraf-nonius.com](http://www.enraf-nonius.com) | [info@enraf-nonius.nl](mailto:info@enraf-nonius.nl) | + 31 - ( 0 ) 10 20 30 600



PARTNER FOR LIFE