

The materials and information provided via this brochure are distributed internationally, but not all the products or services found on this brochure are necessarily available in your country or area. Please contact local office or agents to make sure of the availability of the relevant products or services.

Mindray Medical Netherlands B.V.

Drs.W.van Royenstraat 8, 3871 AN Hoevelaken, The Netherlands
Voice: (31-33) 254-4911
Fax: (31-33) 253-4280

Mindray Medical Italy S.r.l.

Via Leonardo DA Vinci, 158, 20090 Trezzano Sul Naviglio (MI), Italy
Voice: (39-02) 5737-401
Fax: (39-02) 5523-2018

Mindray Medical España S.L.

Avenida Manoteras 38, Edificio B, Bajo B008. 28050-Madrid Spain
Voice: (34-91) 392-3754
Fax: (34-91) 392-3756

Mindray Medical Germany GmbH

Goebelstrasse 21, 3. OG, 64293 Darmstadt
Voice: (49 6151) 3910 0
Fax: (49 6151) 3910 300

Mindray (UK) Limited

Mindray House, Kingfisher Way, Hinchingsbrooke Business Park,
Huntingdon, Cambs, PE29 6FN
Voice: (44) 1480-416-840
Fax: (44) 1480-436-588

Mindray Medical France SARL

Europarc Creteil 1 allée des cerisiers 94000 CRETEIL, France
Voice: (33-1) 4513-9150
Fax: (33-1) 4513-9151

Mindray Medical Sweden AB

Rissneleden 136, 174 57 Sundbyberg, Sweden
Voice: (46) 8 555 54 121

Mindray Medical Rus Co.Ltd

2 Zvenigorodskaya Street, 13, bld 41, 123022 Moscow, Russia
Voice: (7-499) 553-6036
Fax: (7-499) 553-6039

Mindray Medical Colombia S.A.S

Room No.702, Torre UNIKA, Carrera 9 No.77-67, Bogota, Colombia
Voice: (57-1) 3130892 / 3210916

Mindray do Brasil – Comércio e Distribuição de Equipamentos Médicos Ltda.

Rua Tavares Bastos, 329, Perdizes, São Paulo - SP,
CEP 05012-020
Voice: (55-11) 3124-8026
Fax: (55-11) 3078-8035

Mindray Medical Mexico S. de R.L. de C.V.

Félix Parra # 175 Colonia. San José Insurgentes Sur Delegación Benito Juárez 03900 México, D. F.
Voice: (52-55) 5661-9450 / (52-55) 5662-6620
Fax: (52-55) 5662-6597

PT.Mindray Medical Indonesia

The East 11th Floor, Unit 08 Jl. Dr. Ide Anak Agung Gde Agung Kav.
E 3.2 No.1 Jakarta 12950
Voice: (6221) 2902-7280
Fax: (6221) 2902-7283

Mindray Medical Egypt

Ground Floor -11 El Imam Metwally El Sharawi Street,
Sheraton Heliopolis - El Nozha, Cairo, Egypt
Voice: (2) 02 22691047
Fax: (2) 02 22691046

Mindray Medical Australia

ADD: 19 Harker St, Burwood VIC 3125, Australia
Voice: (61-3) 9830 6610
Fax: (61-2) 9550 2841

Mindray Medical Nigeria

ADD: 1004 Estate,Victoria Island, Lagos, Nigeria
Voice: (+234) 7030682053

Mindray Medical India Pvt. Ltd. (Delhi Office)

16th floor, Building No. 9-B, DLF Cyber City, DLF Phase- 3,
Gurgaon- 122002, Haryana
Voice: (91-124) 4632-488
Fax: (91-124) 4632499

Mindray Medical India Pvt. Ltd. (Mumbai office)

B-404, City Point, Andheri Kurla Road, Andheri East, Mumbai,
400 059, India
Voice: (91-22) 4020-0000
Fax: (91-22) 4020-0011

Mindray Medikal Teknoloji Istanbul Ltd. Sti.

KÖYALTı MEVKİİ 29 EKİM CAD.N.1 KUYUMCUKENT 2.NO.PLZ.3.KT.
D.1- 9 34520 BAHÇELİEVLER/YENİBOSNA, İSTANBUL, TÜRKİYE
Voice: (90-212) 482 08 77
Fax: (90-212) 482 08 78

Mindray Medikal Teknoloji Istanbul Ltd. Sti. Branch Office Ankara

TURAN GÜNEŞ BULVARI ALEKSANDER DÜBÇEK CAD. NO:18/A
ÇANKAYA/ANKARA, TÜRKİYE
Voice: (90-312) 441 96 92
Fax: (90-312) 441 96 93

Mindray Medical Vietnam Co.,Ltd.

Unit 2, 4th Floor, "B" Building, Broadway Office Park,
102 Nguyen Luong Bang Blvd., Tan Phu Ward, Dist. 7, HCMC, Vietnam
Voice: (84) 8 5413-6970
Fax: (84) 8 5413-6971

Mindray Medical Thailand Limited

1768 Thai summit Tower Room 407, 4th floor, New patchaburee Road,
Kwaeng Bangkokpi, Khet Huai Khwang, Bangkok 10320, Thailand
Voice: (66) 2652-7288
Fax: (66) 2652-7266

Mindray Medical Korea

KNi Tower 9F 823-19 Yeoksam dong, Seoul-ganum, Korea
Tel: 008225688040
Fax: 008225688042



Accessories for operating tables

| | |
|---|--|
| <p>Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China Tel: +86 755 8188 8998 Fax: +86 755 26582680 E-mail: intl-market@mindray.com www.mindray.com</p> | <p>mindray healthcare within reach are registered trademarks or trademarks owned by Shenzhen Mindray Bio-medical Electronics Co., LTD. ©2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice. P/N:ENG-Accessories for operating tables-210210x44P-20160818</p> <div data-bbox="1155 1442 1371 1498"></div> |
|---|--|



Accessories for Operating Tables

Accessories for operating tables from Mindray are characterised by these functionality, flexibility and high reliability. Developed in close collaboration with end users, the accessories fulfill all kinds of surgery positioning demands. They place the patients in an optimum way that offer benefits both for surgeons and patients.

We use high quality stainless steel materials in our operating tables and accessories. In order to get the surgeons access to C arm, all the pads or plates are radiolucent for X-ray check.

To avoid discomfort in sensitive areas of the patients, we provide gel material pads for patient positioning which are similar as human tissue.



Lateral position



Prone position



Lithotomy position



Supine position



01-11

01

General Accessories

15-20

15

Leg Accessories

23-26

23

Accessories for
Gynecology & Urology

31-34

31

Accessories for controlling
and easier application

12

12-14

Head Accessories

21

21-22

Back Accessories

27

27-30

Orthopedic Accessories

35

35-38

Comfort solution
Accessories




General Accessories



0092-30-102646



0092-30-102634 115-016721-00



115-030153-00

Light arm board with clamp

Accessory Name:
Arm board with rotation function in horizontal position, one piece, with integrated clamp, pad & two pieces of fasten belt

Application:
To support and fasten the arm of the patient for avoiding unexpected moving

Compatible with:

- 8500 / 8300 HyBase
- 6100 HyBase
- 3000 HyBase
- 30 UniBase

Heavy arm board (MFP/PU) w/Clamp

Accessory Name:
Arm board with up & down and rotation function, one piece with one heavy radial universal clamp, pad & two pieces of fasten belt (SFC/PU)

Application:
To support and fasten the arm of the patient for avoiding unexpected moving

Compatible with:

- 8500 / 8300 HyBase
- 6100 HyBase
- 3000 HyBase
- 30 UniBase

Carbon fiber arm board, a pair

Accessory Name:
Carbon fiber arm board, a pair

Application:
Carbon fiber arm board, includes a memorised foam pad. Can allow the 1180 cardiac operating plate to provide support for a patient's upper extremity in a supine position

Compatible with:

- 8500 / 8300 HyBase

General Accessories



115-005693-00

Double arm board with clamp

Accessory Name:
Double layer arm support with one piece of heavy radial clamp for lateral position, with the pad

Application:
To support and fasten both of the arms of the patient as a lateral position to avoid unexpected movement

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



0092-30-102657

Light quadrate clamp

Accessory Name:
Universal clamp (Light, one-piece, quadrate)

Application:
To attach to the side rail for connecting accessories with the table

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



0092-30-102648

Arm board for lateral position w/ clamp

Accessory Name:
Lateral position arm support with onepiece of radial universal clamp, with the pad, one-piece

Application:
To support and fasten the upper arm of the patient in lateral position

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



0092-30-102655

Heavy quadrate clamp

Accessory Name:
Universal clamp (Heavy, one-piece, quadrate)

Application:
To attach to the side rail for connecting accessories with the table

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



115-007047-00

Arm operating table

Accessory Name:
Arm operating table

Application:
To support the arm of the pantient and the hand of the surgeon

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



0092-30-102658

Light radial clamp




Accessory Name:
Universal clamp (Light, one-piece, radial)

Application:
To attach to the side rail for connecting accessories with the table

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase

General Accessories

| | | |
|---|---|--|
|  | <p>Heavy radial clamp</p> <p>Accessory Name: Universal clamp (Heavy, one-piece, radial)</p> <p>Application: To attach to the side rail for connecting accessories with the table</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div><div>6100 HyBase</div><div>3000 HyBase</div><div>30 UniBase</div></div> |
|  | <p>Light Ana. frame with clamp</p> <p>Accessory Name: Anesthesia Frame, one-piece with light quadrate clamp</p> <p>Application: To hang the anesthesia screen during the surgery</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div><div>6100 HyBase</div><div>3000 HyBase</div><div>30 UniBase</div></div> |
|  | <p>Heavy Ana. frame with clamp</p> <p>Accessory Name: Anesthesia Frame , "L" shape design with two adjustable tubes, one-piece with clamp</p> <p>Application: To hang the anesthesia screen during the surgery</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div><div>6100 HyBase</div><div>3000 HyBase</div><div>30 UniBase</div></div> |

0092-30-102656

0092-30-102645

0092-30-102654

| | | |
|---|---|--|
|  | <p>U-Shaped Ana. frame w/ clamps</p> <p>Accessory Name: "U" shape Anesthesia Frame with clamps</p> <p>Application: To hang the anesthesia screen during the surgery</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div><div>6100 HyBase</div><div>3000 HyBase</div></div> |
|  | <p>Suspended arm strap</p> <p>Accessory Name: Suspended arm strap</p> <p>Application: To attach to anesthesia frame for hanging the hand of the patient</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div><div>6100 HyBase</div><div>3000 HyBase</div><div>30 UniBase</div></div> |
|  | <p>Wristlet with clamp</p> <p>Accessory Name: Wristlet with one-piece of quadrate universal clamp</p> <p>Application: To attach to the side rail for fastening the hand of the patient</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div><div>6100 HyBase</div><div>3000 HyBase</div><div>30 UniBase</div></div> |

115-005689-00

115-007051-00

115-007053-00

General Accessories



0092-30-102669 0092-30-102638

Heavy body support with clamp, pubis/sacrum

Accessory Name:
Heavy body support with clamp, pubis
Heavy body support with clamp, sacrum

Application:
To support the body of the patient for avoiding unexpected movement

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



0092-30-102637 0092-30-102671

Heavy shoulder support with clamp, left/right

Accessory Name:
Left/right shoulder support with pad (one piece, Heavy version with one piece of heavy quadrate universal clamp and pad)

Application:
To support the shoulder of the patient for avoiding unexpected movement

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



0092-30-102650

Light body support w/ clamps, a pair

Accessory Name:
Two pieces of body support with a pair of quadrate universal clamps and the pad

Application:
To support the body of the patient for avoiding unexpected movement

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



0092-30-102649

Light shoulder support w/ clamps, a pair

Accessory Name:
Shoulder support (left and right, pair with quadrate universal clamps and pad)

Application:
To support the shoulder of the patient for avoiding unexpected movement

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



115-038996-00

Arm posturing device

Accessory Name:
Arm posturing device

Application:
During surgeries, it is used to allow support for patient's arm to effectively protect patient's brachial nerves

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase



0092-30-102640

Body strap

Accessory Name:
Body strap, one-piece

Application:
To fasten the body of the patient to the table during the surgery

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase

General Accessories



0092-30-102639

Pad for intervertebral disc operations

Accessory Name:
Intervertebral disc operation
pad (Memorise foam, one piece)

Application:
To support the body in
intervertebral operating position

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



115-033429-00

Side rail extension

Accessory Name:
Side rail extension

Application:
Can extend the side rail toward the two ends of the
operating table, to provide more space for installation
of other accessories, such as leg support and foot
plate

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



115-039082-00

Side rail extension

Accessory Name:
Side rail extension (HySRE-2)

Application:
Can extend the side rail towards the two sides of the
operating table, to provide more space for installation
of other accessories

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



115-038243-00

Strap clamps, a pair

Accessory Name:
Strap clamps, a pair

Application:
Can fix the body strap at a specified location on the
side rail

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



0092-30-102641

Cassette fixer with clamp

Accessory Name:
X-ray cassette fixer (compatible with multi-size)

Application:
To fix the X-ray cassette under the table

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



0092-30-102643

Infusion holder with clamp



Accessory Name:
Infusion pole with one piece of heavy radial clamp

Application:
To hang the infusion bag in the procedure

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase




General Accessories

| | | |
|--|--|--|
|  | <p>Table width extender with clamps, a pair</p> <p>Accessory Name: Table width extender with clamps, a pair</p> <p>Application: During surgeries, it is used to widen the operating table top and provide comfortable lying for larger-sized patients</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div><div>6100 HyBase</div><div>3000 HyBase</div><div>30 UniBase</div></div> |
|  | <p>Hand grips with clamps, a pair</p> <p>Accessory Name: Hand grips, a pair</p> <p>Application: Applicable to obstetric operations, for women being in labour to grip.</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div><div>6100 HyBase</div><div>3000 HyBase</div><div>30 UniBase</div></div> |

115-030814-00

115-030812-00

Head Accessories

| | | |
|---|---|--|
|  | <p>Head plate; Head plate, CF</p> <p>Accessory Name: Head plate, single-joint, carbon fibre (optional), with special foam core pad.</p> <p>Application: During surgeries, it is used to provide support to the patient's head</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div></div> |
|  | <p>2-Joint head plate; 2-Joint head plate, CF</p> <p>Accessory Name: Head plate, dual-joint, (optional: carbon fibre), with special foam core pad</p> <p>Application: During surgeries, it is used to provide support to the patient's head with a double-joint design</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div></div> |
|  | <p>Ophthalmic operating plate</p> <p>Accessory Name: Ophthalmic operating head plate (compatible with the adaptor of head position)</p> <p>Application: To support the head of the patient and support the arm of the surgeon during the surgery</p> | <p>Compatible with:</p> <div><div>8500/8300 HyBase</div><div>6100 HyBase</div><div>3000 HyBase</div><div>30 UniBase</div></div> |

115-028574-00 115-028576-00

115-028705-00 115-028703-00

115-007045-00

Head Accessories



0092-30-102642

Adaptor for head position

Accessory Name:
Neurosurgical head support adapter

Application:
To connect the adjustable horseshoe shaped head support (115-005695-00), ophthalmic operating plate (115-007045-00), aluminum headrest system (DJ15-0101-001) and radiolucent headrest system (DJ15-0101-002) with the table

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase



115-005695-00

Head rest, Type1

Accessory Name:
Adjustable horseshoe shaped head rest

Application:
To support the head of the patient providing a stable table for neurosurgery operation. you should select the neurosurgical head support adapter (0092-30-102642) and neurosurgical head support adapter,type 2 (115-034319-00) simultaneously

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



115-012835-00

Head rest, Type 2

Accessory Name:
Adjustable horseshoe-shaped head rest with two installation bars

Application:
To support the head of the patient providing a stable table for Neurosurgery operation. You should select the Neurosurgical head support adapter simultaneously

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



115-034319-00

Universal adaptor for head position

Accessory Name:
Universal adaptor for head position

Application:
To connect the adjustable horseshoe shaped head support(115-005695-00), ophthalmic operating plate(115-007045-00), aluminum headrest system(DJ15-0101-001) and radiolucent headrest system(DJ15-0101-002) with the table

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



Leg Accessories



115-033430-00

Leg restraint cuff, thigh, a pair

Accessory Name:

Leg restraint cuff, thigh, a pair

Application:

Used to secure the patient's thigh.

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



115-033431-00

Leg restraint cuff, shank, a pair

Accessory Name:

Leg restraint cuff, shank, a pair

Application:

Used to secure the patient's shank.

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



115-035792-00 115-037739-00

One-section leg plate, hook ver; One-section leg plate, hook, CF

Accessory Name:

One-section leg plate, hook ver
One-section leg plate, hook, CF




Application:




During surgeries, it is an one-section plate used to support the patient's legs.

Compatible with:

8500/8300
HyBase

Leg Accessories

| | | |
|---|---|--|
|  | 1-joint leg plate, a pair ; 1-joint leg plate, a pair, CF | |
| | <p>Accessory Name:</p> <p>Leg plate, single-joint, (carbon fibre), with special foam core pad, a pair</p> <p>Application:</p> <p>During surgeries, it is used to provide support to the patient's leg</p> | <p>Compatible with:</p> <div>8500 / 8300 HyBase</div> |
| 115-028580-00 115-028578-00 | | |
|  | Rectal positioning device | |
| | <p>Accessory Name:</p> <p>Rectal positioning device</p> <p>Application:</p> <p>To support the coxa and knees of the patient in rcetal surgery</p> | <p>Compatible with:</p> <div>6100 HyBase</div> <div>3000 HyBase</div> |
| 115-007049-00 | | |
|  | 2-joint leg plate, a pair ; 2-joint leg plate, a pair, CF | |
| | <p>Accessory Name:</p> <p>Leg plate, dual-joint abduction, with special foam core pad, a pair (optional: carbon fibre).</p> <p>Application:</p> <p>To support the legs of the patient during the surgery. During surgeries, it is used to provide support to the patient's leg. A double-joint design, helps doctors get closer to the surgical site.</p> | <p>Compatible with:</p> <div>8500 / 8300 HyBase</div> |
| 115-028582-00 115-028584-00 | | |

| | | |
|---|--|--|
|  | One-section leg plate, PU/memory foam pad | |
| | <p>Accessory Name:</p> <p>One section leg plate for PU or SFC (Special Form Core) mattress</p> <p>Application:</p> <p>To support the leg with one-section plate during the surgery</p> | <p>Compatible with:</p> <div>6100 HyBase</div> <div>3000 HyBase</div> <div>30 UniBase</div> |
| 115-005683-00 115-005685-00 | | |
|  | Thigh holder with clamps, a pair | |
| | <p>Accessory Name:</p> <p>Fixation for femur, a pair, with two radial universal clamps</p> <p>Application:</p> <p>To support the femur of the patient during the surgery</p> | <p>Compatible with:</p> <div>8500 / 8300 HyBase</div> <div>6100 HyBase</div> <div>3000 HyBase</div> <div>30 UniBase</div> |
| 115-012841-00 | | |
|  | Shank holder with clamps, a pair | |
| | <p>Accessory Name:</p> <p>Fixation for lower leg, a pair, with two radial universal clamps</p> <p>Application:</p> <p>To support the lower leg of the patient during the surgery</p> | <p>Compatible with:</p> <div>8500 / 8300 HyBase</div> <div>6100 HyBase</div> <div>3000 HyBase</div> <div>30 UniBase</div> |
| 115-012843-00 | | |

Leg Accessories



0092-30-102636 0092-30-102667

Foot plate with clamp, left/right

Accessory Name:

Left/right foot support with one piece of heavy quadrate universal clamp and pad

Application:

To support the foot of the patient for avoiding unexpected movement

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



115-033425-00

Transfer board, hook version

Accessory Name:

Transfer board, hook version lithotomy position.

Application:

During surgeries, it is used to temporarily support the patient's lower extremities when positioning the patient or transferring the patient, and to optimize the pre-operation and post-operation workflows at lithotomy position.

Compatible with:

8500/8300
HyBase



115-033416-00

Meniscus positioning device with clamp

Accessory Name:

Meniscus positioning device with clamp

Application:

To provide support for patient's thigh or foot during knee joint surgeries

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



Back Accessories



115-033428-00 115-038564-00

Short upper back plate ; Short upper back plate, CF

Accessory Name:

Short upper back plate
Short upper back plate , CF

Application:

During Surgeries, provide support to the patient's back

Compatible with:



115-029430-00 115-029432-00

Extension plate,w/ MFP,CF

Accessory Name:

Extension plate, (optional: carbon fibre), with gynaecological cutout, with special foam core pad, length of 250mm.

Application:

With gynecological cutout, can extend the head end or the foot end of the operating table top up to 250mm, to provide more perspective space or more standing space for doctors. Applicable to trauma, spine, gynecology, or urology.
Includes a optional memorised foam pad.

Compatible with:



115-028586-00 115-028588-00

Upper back plate; Upper back plate, CF


Accessory Name:

Upper back plate, (carbon fibre), with special foam core pad.

Application:

During surgeries, provides support to the patient's back.

Compatible with:



115-029434-00

Cardiac operating plate,CF

Accessory Name:

Cardiac operating plate, carbon fibre, with special foam core pad, length of 1,180mm.

Application:

Cardiac operating plate, can provide 360° X-ray radiolucency, extend the foot end of the table top for 1180mm, allows for cardiovascular surgeries and spine surgeries. Includes a memorised foam pad and stand bars.

Compatible with:



115-030151-00

Accessory adaptor,CF

Accessory Name:

Accessory adaptor for Cardiac Operating Plate, carbon fibre, with side rails.

Application:

Provide an additional side rail for the 1180 cardiac operating plate (EU Standard), can be used to install accessories such as arm boards and shoulder supports.

Compatible with:



115-033414-00

Back plate for shoulder operation

Accessory Name:

Back plate for shoulder operation, with special foam core pad.

Application:

During surgeries, it is used to support the patient's back and head; the left and right shoulder plates are detachable, and the installing location at the head can be adjusted upon the patient's height and the shoulder size.

Compatible with:





Accessories for Gynecology & Urology



115-022794-00

Stirrups with clamps, a pair

Accessory Name:

Leg support (Stirrup type), a pair, with clamp

Application:

To support and fasten the leg of the patient

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



115-005697-00

Drainage bowl, HyBase 6100

Accessory Name:

Drainage bowl for Gynecology & Urology

Application:

To drain the liquid in lithotomy position

Compatible with:

6100
HyBase



115-013018-00

Drainage bowl, HyBase 3000

Accessory Name:

Drainage bowl for Gynecology & Urology

Application:

To drain the liquid in lithotomy position

Compatible with:

3000
HyBase

Accessories for Gynecology & Urology



0092-30-102635

Heavy leg support with clamp

Accessory Name:

One piece of heavy version leg support, with one piece of heavy radial universal clamp, pad & one piece of fasten belt

Application:

To support and fasten the leg of the patient

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



0092-30-102647

Light leg support with clamps, a pair

Accessory Name:

Light leg support, a pair, with radial universal clamp

Application:

To support and fasten the leg of the patient

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



115-007043-00

Leg holder with loops, with clamps, a pair

Accessory Name:

A pair of Leg holders with loops, with radial universal clamps

Application:

To support and fasten the leg of the patient

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



115-034318-00

Drainage bowl with clamp, universal

Accessory Name:

Universal drainage bowl

Application:

To drain the liquid in lithotomy position

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase



115-030812-00

Hand grips with clamps, a pair

Accessory Name:

Hand grips, a pair

Application:

Applicable to obstetric operations, for women being in labor to grip.

Compatible with:

8500/8300
HyBase

6100
HyBase

3000
HyBase

30
UniBase

Orthopedic Accessories



115-029436-00

Ortho extension device w/ cart,CF

Accessory Name:

Ortho extension device, carbon fibre, including the cart

Application:

During the orthopedic traction surgery, it is used to pull and fix the affected limb to meet the need of x-rays; the device is made of highly-transmissible carbon fibre materials to improve image qualities.

Compatible with:



115-007055-00

Tibia & fibula extension assembly

Accessory Name:

Tibia extension assembly

Application:

To provide tibial and fibula traction

Compatible with:



115-033422-00

Upper arm extension device with clamps

Accessory Name:

Upper arm extension device

Application:

During surgeries, it is used to fix and pull patient's forearms; consists of auxillary pillar, elbow joint pillar, palm & finger fixer and tractor

Compatible with:



115-033409-00

Ortho extension w/ cart, hook version

Accessory Name:

Ortho extension device, hook version, stainless steel, including the cart

Application:

To provide leg traction during orthopedic surgery

Compatible with:



115-033421-00

Forearm extension device with clamps

Accessory Name:

Forearm arm extension device

Application:

During surgeries, it is used to fix and pull patient's forearms; consists of auxillary pillar, elbow joint pillar, palm & finger fixer and tractor

Compatible with:



115-033410-00

Ortho extension w/ cart, UniBase 30

Accessory Name:

Ortho extension w/ cart, UniBase

Application:

To provide leg traction during orthopedic surgery

Compatible with:



Orthopedic Accessories



115-009425-00

Orthopedic extension device with cart

Accessory Name:

Orthopedic extension device, including the cart

Application:

To provide leg traction during orthopedic surgery

Compatible with:

- 6100
HyBase
- 3000
HyBase



115-011817-00

Pelvis width extender, a pair

Accessory Name:

Triangle support board, 2 pieces.

Application:

During the orthopedic traction surgery, it is used to widen the layer board to provide comfortable support.

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase



115-007057-00

Lateral femur extension assembly

Accessory Name:

Lateral femoral extension assembly

Application:

During the orthopedic traction surgery, it is used to pull and fix the affected limb at supine position to meet the need of x-rays.

Adaptive accessories:

Light arm board, heavy arm board, arm board for lateral position, double arm board, heavy radial clamps.

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase



115-033415-00

Upper arm operating table

Accessory Name:

Upper arm operating table

Application:

To support the upper arm of the patient.

Compatible with:

- 8500/8300
HyBase
- 6100
HyBase
- 3000
HyBase
- 30
UniBase



115-033413-00

Arthroscopic leg holder with clamp

Accessory Name:

Arthroscopic leg holder with clamp

Application:

To fasten and support the knee of the patient

Compatible with:

- 6100
HyBase
- 3000
HyBase
- 30
UniBase

Accessories for controlling and easier application



115-012837-00

Mayo instrument tray with clamp

Accessory Name:

Mayo instrument tray

Application:

To place surgical instrument.

Compatible with:

- 8500/8300 HyBase
- 6100 HyBase
- 3000 HyBase
- 30 UniBase



Comfort solution Accessories

Unit: cm



- 115-005822-00 Ankle protection strap, 45×7×1



- 115-005826-00 Table pads



- 115-005900-00 Universal head pad, large, 25×23×6
- 115-005832-00 Universal head pad, small, 25×18×3



- 115-005864-00 Prone position pillow, large, 28×24×14
- 115-005830-00 Prone position pillow, small, 25×20×10
- 115-038074-00 Prone pillow, full support, 23×20×11
- 115-038075-00 Paediatric prone pillow, full 21×18×11



- 115-005872-00 Head ring, large, 20×7×7
- 115-005873-00 Head ring, middle, 20×7×5
- 115-005874-00 Paediatric head ring, large, 14×5×3
- 115-005875-00 Paediatric head ring, small, 10×4×2.5
- 115-005876-00 Infant head ring, 8×3×2
- 115-005877-00 Head ring, small, 20×7×3
- 115-005878-00 Paediatric head ring, middle, 14×7×2.5



- 115-005885-00 Horseshoe head ring, large, 20×7×7
- 115-005886-00 Horseshoe head ring, middle, 20×7×5
- 115-005887-00 Paediatric horseshoe head ring, L14×5×3
- 115-005888-00 Paediatric horseshoe head ring, S10×4×2.5
- 115-005889-00 Infant head ring, 8×3×2
- 115-005890-00 Horseshoe head ring, small, 20×7×3

Unit: cm



- 115-005892-00 Universal supine pad, 18×18×8



- 115-005896-00 Ophthalmic head pad, large, 23×23×4
- 115-005897-00 Ophthalmic head pad, small, 22×19×6
- 115-005898-00 Pillow for supine position, 24×22×7
- 115-038076-00 Ophthalmic pad, paediatric, 15×15×3.5



- 115-005818-00 Universal square pad-small, 30×20×1.5
- 115-005819-00 Universal square pad-middle, 40×24×1.5
- 115-005820-00 Universal square pad-large, 50×30×1.5
- 115-038078-00 Universal hip pad, 50×40×1.5



- 115-005903-00 Universal arm board pad, short, 40×16×1
- 115-005904-00 Universal arm board pad, long, 60×16×1



- 115-005906-00 Elbow protection strap, 40×15×1.5



- 115-005908-00 Wrist protector, 25×5×1

Comfort solution Accessories

Unit: cm



- 115-005910-00 Lateral positioner, 72×46×15



- 115-005788-00 Pad for intervertebral disc, 63×46×15
- 115-038080-00 Prone position pad, hard inner, 60×49×18cm



- 115-038077-00 Struma pad, 48×39×9



- 115-005794-00 Pad for contertraction post, 27×20×1.5



- 115-005803-00 Short gel body roll, large, 30×15×7
- 115-005804-00 Short gel body roll, small, 30×10×4
- 115-005805-00 Long gel body roll, large, 50×15×7
- 115-005806-00 Long gel body roll, small, 50×10×4
- 115-005807-00 Short gel body roll, middle, 30×10×7
- 115-005808-00 Long gel body roll, middle, 50×10×7
- 115-038079-00 Gel body roll, paediatric, 20×10×7



- 115-005796-00 Body strap, 60×7×1

Unit: cm



- 115-005811-00 Arm/leg pad, long, 50×15×4
- 115-005812-00 Arm/leg pad, short, 40×15×4



- 115-005814-00 Heel pad, one piece, 18×10×7



- 115-005790-00 Chest-hipbone pad, 50×20×7
-
-

A close-up photograph of a dark, textured carbon fiber table top. A polished metal hinge is attached to the edge. The image is part of a larger graphic with a red background on the right.

Carbon Fiber table top

mindray
healthcare within reach