

OTOPLASTIK
OTOPLASTY
OTOPLASTIA
OTOPLASTIE
OTOPLASTICA

SET

ZS 205/01

Siebvorschlag - Suggested instrument set



ZS 205/01	Content:	Qty
AB 060/16	STANDARD TISSUE FORCEPS, 1X2 T., STRAIGHT, 16,0 CM	1
AB 116/09	ADSON-BROWN TISSUE FORCEPS, 9X9 T., STRAIGHT, 12,0 CM	1
AC 981/13	GOLDMAN-FOX SCISSORS TUC, CURVED, 13,0 CM	1
AD 010/03	SCALPEL HANDLE NO.3	2
AD 075/10	PCK=100 SCALPEL BLADES, FIG.10	1
AE 466/13	HALSEY NEEDLE HOLDER TUC, 13,0 CM	1
AK 022/15	KILNER DELICATE HOOK FINE, 15,0 CM	1
AK 072/02	JOSEPH DELICATE HOOK, 2 PRONGS FINE, 16,0 CM	1
KO 409/00	COTTLE(MILLARD) DELICATE HOOKS, 2-PRONGS	1
KU 300/11	MIXING BOWL, 116 x 50 MM	1
XR 830/03	CONNECTION CORD BIPOLAR, 300 CM	1
XR 894/16E	BI-COAGULATION-FORCEPS, NEEDLE 16,0 CM	1
KYE 304/50	WIRE BASKET, 267 x 125 x 50 MM	1
KYE 810/01	SILICON MAT, 275 x 125 MM	1


Sterilizing Container KYE 351/01, if required.



STANDARD 1:2
AB 060/16
16,0 cm/6¹/₄"



ADSON-BROWN
AB 116/09
12,0 cm/4³/₄"



GOLDMAN-FOX
AC 981/13
13,0 cm/5"



AD 010/03
No. 3



HALSEY
AE 466/13
13,0 cm/5"



KILNER
AK 022/15
15,0 cm/6"



JOSEPH
AK 072/02
16,0 cm/6¹/₄"

neck lift

rhino basic

rhino advanced

otoplasty

facelift

facelift endoscopy

blepharoplasty

hair transplantation



COTTLE (MILLARD)
KO 409/00



KU 300/11
Ø x h
116 x 50 mm 0,35 l

bipolar
Connecting cable



MARTIN / BERCHTHOLD
XR 830/03
300 cm



XR 894/16E
16,0 cm/6¹/₄"



KYE 304/50

L(mm)	B(mm)	H(mm)
267	125	50



KYE 810/01

L(mm)	B(mm)
275	125