


| Contents: | | Qty |
|-----------|---|-----|
| AA 150/14 | HALSTED-MOSQUITO ART.FCPS STR 14,0 CM | 2 |
| AA 151/14 | HALSTED-MOSQUITO ART.FCPS CVD 14,0 CM | 2 |
| AA 711/09 | JONES TOWEL FORCEPS 9,0 CM | 2 |
| AB 628/09 | ZIEGLER CILIA FORCEPS 9,0 CM | 1 |
| AB 785/01 | GRAEFE FIX.FCPS W/O SPRING CATCH 11,0 CM | 1 |
| AB 816/01 | CASTROVIEJO SUTURE FORCEPS 10,0 CM | 2 |
| AB 824/02 | BONACCOLTO SUTURE FORCEPS 10,0 CM FINE | 1 |
| AB 829/01 | HARMS TYING FORCEPS STRAIGHT 10,0 CM | 1 |
| AB 840/01 | DESMARRES CHALAZION FORCEPS FIG.1 9,0 CM | 1 |
| AB 841/01 | LAMBERT CHALAZION FORCEPS 9,0 CM ROUND | 1 |
| AB 844/20 | BERKE PTOSIS FORCEPS 11,0 CM, 20 MM | 1 |
| AC 450/10 | IRIS SCISSORS STR SH/SH 10,5 CM | 1 |
| AC 475/10 | STEVENS TENOTOM.SCS POINTED CVD 10,5 CM | 1 |
| AC 782/11 | WESTCOTT TENOTOMY SCISSORS BLUNT 11,0 CM | 1 |
| AC 797/01 | WESTCOTT STRABISMUS SCS STR SHARP 11,5 CM | 1 |
| AE 663/14 | CASTROVIEJO NH CVD RA SERRATED 14,0 CM | 1 |
| AE 688/14 | KALT MICRO-NEEDLE HOLDER STR 14,0 CM | 1 |
| AK 074/01 | GUTHRIE DELICATE HOOK FIG.1 16,0 CM | 1 |
| AK 078/13 | KNAPP DELICATE HOOK BLUNT 13,0 CM | 1 |
| AK 090/12 | DESMARRES RETRACTOR 12 MM 13,0 CM | 1 |
| AK 090/14 | DESMARRES RETRACTOR 14 MM 13,0 CM | 1 |
| AK 090/16 | DESMARRES RETRACTOR 16 MM 13,0 CM | 1 |
| AR 282/20 | CASTROVIEJO CIRCLE 8,0 CM STRAIGHT | 1 |
| KP 015/01 | GRAEFE STRABISMUS HOOK #1 13,0 CM | 1 |
| KP 015/02 | GRAEFE STRABISMUS HOOK #2 14,0 CM | 1 |
| KP 248/11 | JAEGER LID PLATE 11,0 CM | 1 |
| KP 325/05 | BONN SCLERA KNIFE 5MM 12,5 CM | 1 |
| KP 415/10 | MEYHOEFER CHALAZION CURETTE 13,0 CM | 1 |
| KP 415/20 | MEYHOEFER CHALAZIONCURETTE 13,0 CM 2,0 MM | 1 |
| KP 415/25 | MEYHOEFER CHALAZIONCURETTE 13,0 CM 2,5 MM | 1 |
| KP 415/35 | MEYHOEFER CHALAZIONCURETTE 13,0 CM 3,5 MM | 1 |
| KP 700/00 | BISHOP-HARMON ANTERIOR CHAMBER | 1 |
| SK 305/05 | WIRE BASKET, 485 X 255 X 50 MM | 1 |

Sterilization Container KYE 200/61, if required



AB 785/01
ohne Schloss
without spring catch
sin cierre de resorte
sans arrêt à ressort
senza arresto a molla

HALSTED-MOSQUITO
AA 150/14 AA 151/14
14,0 cm/5 1/2"

JONES
AA 711/09
9,0 cm/3 1/2"

ZIEGLER
AB 628/09
9,0 cm/3 1/2"

GRAEFE
11,0 cm/4 1/3"

CASTROVIEJO 1:2
AB 816/01
0,12 mm
10,0 cm/4"

BONACCOLTO
AB 824/02
10,0 cm/4"



HARMS
AB 829/01
10,0 cm/4"



DESMARRES
AB 840/01
9,0 cm/3 1/2"




LAMPERT
AB 841/01
9,0 cm/3 1/2"




BERKE
AB 844/20
11,0 cm/4 1/3"
20 mm

Irisscheren
Iris Scissors
Tijeras para iridectomía
Ciseaux à iridectomie
Forbici per iridectomia



AC 450/10
10,5 cm/4¹/₄"

Sehnenscheren
Tenotomy Scissors
Tijeras para tenotomía
Ciseaux à ténatomie
Forbici per tenotomia



STEVENS
AC 475/10
10,5 cm/4¹/₄"



WESTCOTT
AC 782/11
11,0 cm/4¹/₃"



WESTCOTT
AC 797/01
11,5 cm/4¹/₂"



CASTROVIEJO
AE 663/14
14,0 cm/5¹/₂"



KALT
AE 688/14
14,0 cm/5¹/₂"

spitz
pointed
agudas
pointues
acuta

gebogen
curved
curvas
courbes
curve

Messbereich 20 mm
measuring range 20 mm
campo de medición 20 mm
étendue de mesurage 20 mm
ampiezza di misura 20 mm



GUTHRIE
AK 074/01
16,0 cm/6¹/₄"



KNAPP
AK 078/13
13,0 cm/5"



DESMARRES
13,0 cm/5"



AK 090/12
12 mm



AK 090/14
14 mm



AK 090/16
16 mm



CASTROVIEJO
AR 282/20
8,0 cm/3"




GRAEFE
KP 015/01
13,0 cm/5"
Fig. 1



GRAEFE
KP 015/02
14,0 cm/5¹/₂"
Fig. 2




JAEGER
KP 248/11
11,0 cm/4 1/2"



BONN
KP 325/05
12,5 cm/ 5"
5 mm




MEYHOEFER
13,0 cm/5"



KP 415/10
Ø 1,0 mm

KP 415/20
Ø 2,0 mm



KP 415/25
Ø 2,5 mm

KP 415/35
Ø 3,5 mm

Set komplett
Set complete
Juego completo
Jeu complet
Set completo



BISHOP-HARMON
KP 700/00



| L(mm) | B(mm) | H(mm) | |
|-------|-------|-------|-----------|
| 485 | 255 | 50 | SK 305/05 |