


Contents:		Qty
AA 060/03	DIEFFENBACH BULLDOG CLAMPS STR 3,5 CM	2
AA 114/10	HARTMANN FORCEPS STR 10,0 CM	2
AA 150/12	HALSTED-MOSQUITO FCPS STR 12,5 CM	2
AA 711/09	JONES TOWEL FORCEPS 9,0 CM	2
AA 751/13	BACKHAUS TOWEL FORCEPS 13,0 CM	2
AB 710/10	GRAEFE IRIS FORCEPS STR SERR. 10,0 CM	1
AB 750/07	BARRAQUER-KATZIN FIX.FCPS 1X2 T. 7,0 CM	1
AB 816/01	CASTROVIEJO SUTURE FORCEPS 10,0 CM	2
AB 824/02	BONACCOLTO SUTURE FORCEPS 10,0 CM FINE	1
AB 825/02	BONN SUTURE FORCEPS 10,0 CM ULTRA FINE	2
AB 848/11	FAULKNER LENS HOLDING FORCEPS 11,0 CM	1
AC 466/10	IRIS SCISSORS STR W LARGE BOWS 10,0 CM	1
AC 467/10	IRIS SCISSORS CVD W LARGE BOWS 10,0 CM	1
AC 474/11	STEVENS TENOTOM.SCS.STR POINTED 11,5 CM	1
AC 776/01	VANNAS MOD.TÜBINGEN STRAIGHT 8,5 CM	1
AC 780/11	WESTCOTT TENOTOMY SCS.POINTED 11,0 CM	1
AE 657/14	CASTROVIEJO NH CVD SMOOTH 14,0 CM	1
AE 660/14	CASTROVIEJO NH STR RA SMOOTH 14,0 CM	1
AR 282/20	CASTROVIEJO CIRCLE 8,0 CM STRAIGHT	1
KP 002/12	JAMESON MUSCLE HOOK 12,0 CM	2
KP 010/14	GASS RETINAL DETACHMENT HOOK 14,0 CM	2
KP 211/08	LIEBERMAN EYE SPECULA 15 MM 8,0 CM ASP.	1
KP 261/15	SAUTTER ORBITAL SPATUL MEDIUM 15,0 CM	1
KP 262/13	ARRUGA PROTECTOR 13,0 CM	1
KP 354/12	FOREIGN BODY NEEDLE STRAIGHT 12,0 CM	1
KP 390/00	VITRECTOMY SCISSORS STRAIGHT SMOOTH 0°	1
KP 394/00	VITRECTOMY GRASPING FCPS STRAIGHT SMOOTH	1
KP 700/00	BISHOP-HARMON ANTERIOR CHAMBER	1
KP 760/21	INJECTION CANNULA 21G 4,5 CM	1
SK 305/05	WIRE BASKET, 485 X 255 X 50 MM	1

Sterilization Container KYE 200/61, if required



DIEFFENBACH
 AA 060/03
 3,5 cm/1 1/2"




HARTMANN
 AA 114/10
 10,0 cm/4"




HALSTED-MOSQUITO
 AA 150/12
 12,5 cm/5"



JONES
 AA 711/09
 9,0 cm/3 1/2"



BACKHAUS
 AA 751/13
 13,0 cm/5"



GRAEFE
 AB 710/10
 10,0 cm/4"



BARRAQUER-KATZIN 1:2 (Colibri)
 AB 750/07
 7,0 cm/2 3/4"



CASTROVIEJO 1:2
 AB 816/01
 0,12 mm
 10,0 cm/4"



BONACCOLTO
 AB 824/02
 10,0 cm/4"



MOD. BONN
 AB 825/02
 10,0 cm/4"




FAULKNER
 AB 848/11
 11,0 cm/4 1/3"

ZS 160/85

VITREKTOMIE (PARS PLANA) VITRECTOMY (PARS PLANA) VITRECTOMÍA PARS-PLANA
 VITRECTOMIE (PARS PLANA) VITRECTOMIA PARS-PLANA

mit großen Ringen
 with large bows
 con anillos grandes
 avec grands anneaux
 con anelli grandi



AC 466/10
 10,0 cm/4"

spitz
 pointed
 agudas
 pointues
 acuta



STEVENS
 AC 474/11
 11,5 cm/4 1/2"

gerade
 straight
 rectas
 droits
 rette



VANNAS MOD.
 TÜBINGEN
 AC 776/01
 8,5 cm/3 1/3"

gerade
 straight
 rectas
 droits
 rette



WESTCOTT
 AC 780/11
 11,0 cm/4 1/3"



CASTROVIEJO
 AE 657/14
 14,0 cm/5 1/2"



CASTROVIEJO
 AE 660/14
 14,0 cm/5 1/2"

Messbereich 20 mm
 measuring range 20 mm
 campo de medición 20 mm
 étendue de mesure 20 mm
 ampiezza di misura 20 mm



CASTROVIEJO
 AR 282/20
 8,0 cm/ 3"



JAMESON
 KP 002/12
 12,0 cm/4 3/4"



GASS
 KP 010/14
 14,0 cm/5 1/2"


T
 B



LIEBERMAN
 KP 211/08
 8,0 cm/3 1/4"
 B = 15 mm



SAUTTER
 KP 261/15
 15,0 cm/6"




ARRUGA
 KP 262/13
 13,0 cm/5"



KP 354/12
 12,0 cm/ 4³/₄"



BISHOP-HARMON
 KP 700/00




KP 390/00
 0° KP 394/00

ohne Spülung
 without irrigation
 sin irrigación
 sans irrigation
 senza irrigazione



Set komplett
 Set complete
 Juego completo
 Jeu complet
 Set completo



KP 760/21
 4,5 cm/2"
 21G



L(mm)	B(mm)	H(mm)	
485	255	50	SK 305/05